kit and refills



myQuickmat hoop kit

circumferential matrix system

REF 6600

40	pcs.	Quickmat hoop, preloaded circumferential matrices, assorted (20 x regular, 20 x tabbed)
2	pcs.	myHooper instrument
100	pcs.	myHoopy wedge, interdental plastic wedges, assorted

for more information

myHooper instrument

Contents		REF
1 pc.	myHooper instrument	6602



Quickmat hoop

preloaded circumferential matrices



Contents		Size	REF
	50 pcs. Quickmat hoop, regular	small 5.5 mm	6603
	50 pcs. Quickmat hoop, regular	small 6.5 mm	6604
	50 pcs. Quickmat hoop, regular	large 5.5 mm	6605
	50 pcs. Quickmat hoop, regular	large 6.5 mm	6606

Contents		Size	REF
	45 pcs. Quickmat hoop, tabbed	small 5.5 mm	6607
	45 pcs. Quickmat hoop, tabbed	small 6.5 mm	6608
	45 pcs. Quickmat hoop, tabbed	large 5.5 mm	6609
	45 pcs. Quickmat hoop, tabbed	large 6.5 mm	6610
	45 pcs. Quickmat hoop, tabbed	23 small, 22 large 8.0 mm	6611

myHoopy wedge

interdental plastic wedges

Contents			Size	REF
	100 pcs.	myHoopy wedge	yellow (S)	7602
	100 pcs.	myHoopy wedge	green (M)	7603
	100 pcs.	myHoopy wedge	teal (L)	7604
	100 pcs.	myHoopy wedge	blue (XL)	7605
	100 pcs.	myHoopy wedge, assorted	(30xS, 30xM, 20xL, 20xXL)	7601



MKTG-0169-FL-EN REV 00 06.2026

Polydentia SA
via Cantonale 47
6805 Mezzovico-Vira
Switzerland

phone: +41 91 946 29 48
info@polydentia.ch
polydentia.ch

polydentia
swiss manufacture

SWISS MADE

myQuickmat hoop

new!

circumferential control you can count on



ADVANTAGES

- easy placement similar to a Tofflemire retainer — without its drawbacks
- detachable instrument for fast placement, tightening and adjustment
- anatomically contoured band for leak-free margins and natural restorations
- compact design that is ideal for tight spaces, including second and third molars
- 35-micron hard contoured foil: strong, shape-retentive and passes tight contacts easily
- colour-coded adjustable plastic retainer that stays secure without loosening
- swivelling retainer arm for use in any quadrant

polydentia
swiss manufacture

circumferential matrix system

SWISS MADE

new!

myQuickmat hoop

circumferential matrix system





"The Quickmat hoop matrix is extremely user-friendly: even extensive composite restorations can be completed quickly, with excellent predictability of the final result.
Interproximal contacts and tooth anatomy are preserved, and the time required for the final finish is always minimised."

Dr. Thomas Casanova - Switzerland

Clinical image by courtesy of Dr. Perucchi & Dr. Scognamiglio, Switzerland

The myQuickmat hoop kit is a circumferential matrix system that preserves the familiar placement workflow of a Tofflemire retainer, while eliminating many of its traditional constraints. With easy placement, anatomical contours, a compact design, hard foil, and full swivelling capability, myQuickmat hoop sets a new benchmark for modern circumferential matrix band systems.

The system includes a selection of regular and tabbed preloaded circumferential matrices, an instrument and interdental plastic wedges, which are useful for stabilising the matrices and separate the adjacent tooth.

watch the video



Quickmat hoop

The **Quickmat hoop** preloaded matrix is a contoured circumferential matrix designed to deliver precision and consistent clinical results. Its detachable instrument (myHooper instrument) allows fast, intuitive placement using a Tofflemire-style technique but without leaving a bulky holder in the patient's mouth. This eliminates interference with soft tissues and prevents the matrix from being displaced during the procedure.
The colour-coded adjustable retainer remains securely in place without loosening, while the swivelling arm allows use in any quadrant.

Quickmat hoop preloaded matrices are available regular and tabbed, in two circumferential sizes (small and large) and three height options (5.5, 6.5, 8.0mm).

The Quickmat hoop matrices are available in a tabbed version to help with large multi-surface composite restorations. The featured Control Hoop Tabs are specifically designed to facilitate reliable tight contacts when the matrix band unintentionally pulls away from adjacent teeth and open contacts occur.





The shape at the base of the band nips in and then flares out. It imitates the shape of a regular tooth which, obviously, is not straight up and down.



To maintain firm contact and ensure tight interproximal contacts:
- **option 1**
  Pull on the tab while curing the proximal composite.
- **option 2**
  Use flowable composite to tack the tabs to the adjacent tooth.









myHoopy wedge

The **myHoopy wedge** provides reliable separation, effective tissue compression and a secure matrix seal to help prevent seepage. Its parallel walls and gradually increasing rigidity allow for smooth insertion and firm stabilization once in position.

